

SPAN Encore




- GlideAlign® Retractable Deck Technology
- Low-bed design with low position auto-pause
- ReadyWide® integrated for easy width expansion

science, safety, simplicity

 savaria®

SPAN Encore



The GlideAlign Retractable Deck

Keeps the user vertically aligned, prevents sliding during HOB elevation, and reduces repositioning needs.




Enhance levels of care and comfort for the user, plus ergonomics and injury protection for the caregiver

Packed with features:


- Intelligent Smart Stop[®] low position auto-pause improves low-bed safety and ergonomics
- Advanced Clinical Positioning includes Comfort Chair, T/RT, and Cardiac Care
- Hospital-grade ball bearing swivel casters, provide full mobility at any height. Head-end includes directional lock-out.
- The bed travels virtually straight up and down, while also providing ample underbed clearance for your convenience.
- Dual-position, dual intensity under-bed lighting.
- Trapezoidal slat profile offers superior strength, ventilation, and cleaning ease
- Integrated ReadyWide[®] expanders to 42"
- Foot-end controller with staff lockout
- Baseboard Bumper, Fixed Mattress Retainers and Vascular Foot Ratchet

Customize with:

- Bed extenders to 84" length
- Ergonomic, 3-position side assist rail bed for sturdy ingress/egress support
- Molded, half-head assist rails
- Built-in backup battery
- Headboard and footboard options
- Gooseneck, headboard-mounted pendant holder
- Bed & deck angle indicator kit
- Side carrying transport dollies to move beds through narrow halls and ergonomic maintenance
- USB Charging Kit



Choice of Finish & Profile



SPECIFICATIONS

Height Floor to Mattress Support	7-7/8" - 31" / 200 mm - 787 mm
Deck Angles – Back Elevation	0 – 70°
Deck Angles – Knee Elevation	0 – 30°
Bed Length – Overall Length	82", 86", 90" / 2083 mm, 2184 mm, 2286 mm
Sleep Surface – Retraction Travel	12" / 305 mm
Bed Length – Sleep Surface	76", 80", 84" * / 1930 mm, 2032 mm, 2134 mm *
Bed Width – Sleep Surface	36", 39", 42" / 889 mm, 991 mm, 1067 mm
Bed Width – With Side Rails	40-1/2", 43-1/2", 46-1/2" / 1029 mm, 1105 mm, 1181 mm
Weight of Bed (without accessories)	253 lbs / 115 kg
Safe Working Load	600 lbs / 272 kg
Mobility (at any height)	7-7/8" - 31" / 200 mm - 787 mm
Bed End Panels	Permanent or removable mounting
Ingress Protection Rating	IPX4 (Electrical components IP66)
Electric Motors	Four (4) DC self-lubricating linear actuators

* Requires the 4" Bed Extender Kit

WARRANTY

Electrical	5 years
Frame	15 years
Welds	15 years

ITEM CODE DESCRIPTION

MC9RW	Encore with ReadyWide®
MC9RZ	Encore with built-in backup battery
MC9RWZ	Encore ReadyWide® with built-in backup battery
R9308H	1/2 side molded rails head, pair
R9519	R9519 Split side rail pair
R9958-L/R	Single molded Pivot Assist rail, quick width adjustable, frame mounted
R9959	Pair molded Pivot Assist rail, quick width adjustable, frame mounted
Q6731	Trapeze mounting frame
Q6730	Trapeze Pole, offset
BS79(B)	I.V. Pole (adjustable, board mounted)
Q6785	4" Bed extender kit - factory installed
Q6689	Expander kit (39"/42" quick adjustment)
Q6701	Battery Back-up kit - factory installed
Q6699	Transport Dolly Set
Q6729	Foot end equipment bar
Q6531	Roller bumper
CS82	Board mounted base for I.V. Pole
Q6760	Bed Extender Storage Kit

CERTIFICATIONS

ANSI/AAMI ES60601-1:2005 - Medical electrical equipment
CAN/CSA-C22.2 NO. 60601-1-08 - Medical electrical equipment
CAN/CSA-C22.2 NO. 60601-2-52:11 - Medical Beds
CAN/CSA-C22.2 NO. 60601-1-2-08 - EMC



Get more information

800.268.8671(Span)

savariapatientcare.com



Span-America

| 800.268.8671

| PO Box 5231 Greenville SC 29606

Span-America is a division of Savaria Corporation, a global leader in accessibility. Learn more about Savaria at savaria.com

Control Document #