

FEATURED HOSPITAL BEDS



Full Electric Standard/Low Hospital Bed



- 450 lb. weight capacity frame
- Bariatric size available
- Adjustable height with dual casters
- Easy delivery with split spring design; Independent operation of head, foot, and hi/lo motors; 84" length extension kit available.
- Available for rental, rent to own 9 months and purchase



Halo Safety Ring
Available for Install

Prius Healthcare Descend Ultra-Low Floor Hospital Bed



- 375 lb. weight capacity frame
- Enhanced safety with ultra-low floor design for fall prevention
- Electric height adjustment from 2.7" to 25"
- Customizable backrest, footrest, Trendelenburg functions; Durable steel and wood construction
- Battery backup for emergency lowering



Optional Fall Mat

Transfer Master Supernal 5 Hi-Lo Hospital Bed



- 300 lb. weight capacity
- Full electric bed with head, foot, and Hi-low adjustability resembling a regular bed
- Available in Twin, Full, Queen, or Dual King sizes
- Supernal 5 features Trendelenburg, reverse Trendelenburg, cardiac chair positions, variable massage, and backlit wireless hand control
- Adjustable bed deck range of 12"-24"



Optional Bed Rail
and Bamboo Cover

BROWSE ALL



HOSPITAL BEDS
& MATTRESSES

- ✓ *Trained Bed and Mattress Experts to Assist You*
- ✓ *3,000sq ft. Showroom to Test Try Options*
- ✓ *White Glove Delivery and Setup*
- ✓ *Freedom Mobility Promise with Every Purchase*




📍 Hanover, MD Showroom

FEATURED MATTRESSES



Span America Geo-Mattress UltraMax



- 500 lb weight capacity with an 8-year warranty
- Advanced foam mattress with three patented features, including Shear Transfer Zones for preventing pressure points
- Bi-directional stretch cover, fluid-proof top cover, and over 800 articulating cells for enhanced comfort and durability

Prius VLX Foam Mattress




- 500 lb Weight Capacity
- Tri-layered foam mattress with cool visco elastic top layer
- Zoned comfort foam, heel soft foam, and firm support base
- Full-length firm side rails, urethane 4-way stretch top cover
- Available in standard, bariatric, and custom sizes.

Prius Salute RDX Mattress System



- 500 lb Weight Capacity
- Budget-friendly premium quality alternating air mattress
- Therapeutic pressure with power outage inflation capability
- 3-in-1 alternating air cycle, 8 adjustable comfort levels
- Low friction, low shear cover for maximum protection and patient comfort

Span America CustomCare Convertible Mattress



- 500 lb Weight Capacity (non-powered mode), up to 350 lb (powered modes)
- Lateral rotation and alternating pressure in pump operation
- Self-enclosed vertical air cylinders for pump-free use
- Easily converts to a powered surface with a convenient control unit
- Shear transfer zone technology for skin protection