

AURA™

Premium Home Hospital Bed



SonderCare
Healthy Living Adjustable Beds





“ SonderCare has given my mom her independence back. It’s so much easier for her to get out of bed and it’s really comfortable. She’s so proud of her bed.”
M. Marisek

Bed Positions



BACKREST



KNEES



CHAIR



ZERO GRAVITY



TRENDELENBURG & REVERSE



HIGH-LOW



One-Button Chair Position for Comfortable Eating, Reading or Watching TV



Maximum High Height for Better Care and Repositioning




Illuminating Hand Control for High Visibility Even in Dark Settings



Zero Gravity Position to Help Alleviate Pressure and Back Pain



Multi-Height Assist Rails for Security and Assistance During Transfers



Hidden Casters and One Central Brake Pedal for Stability and Easy Movement



Ultra Low Height for Easy Transfers and Fall Prevention



The Aura™ Bed Size & Mattress Selection

The Aura™ Bed Available in Platinum or Premium



Aura™ Premium 39" with Comfort™ Mattress



Aura™ Premium 39" with Dream™ Mattress



Aura™ Premium Extra Wide 48" with Dream™ Mattress



Aura™ Premium Espresso Woodgrain



Aura™ Platinum Slate Gray Fabric



Aura™ Platinum Shown with two rails and Dream™ Mattress



Aura™ Platinum Shown with no rails and Comfort™ Mattress

The Extra Wide Aura™ models are 48" mattress width, providing users with a little extra space along with all the safety and mobility functions of the original Aura™. Available in both *Premium* and *Platinum* versions.

The Aura™ Platinum is the perfect blend of advanced features and beautiful design. Upholstered in fluid-proof and stain-resistant Crypton™ ultra-soft fabric.

All SonderCare mattresses are designed to work with your bed system to provide optimal comfort and pressure redistribution


Dream™ Mattress



Gel-Infused Visco-Elastic Memory Foam Reduces Pressure and Regulates Temperature




Dream™ Bamboo Quilt-Top Mattress is Reversible for Firm or Plush Sides



Comfort™ Mattress includes Memory Foam without Quilt-Top (Non-Reversible)



All Mattresses Are Outfitted with 100% Fluid-Proof Zippered Stretch Covers



Multi-Zone Design Provides Durability, Comfort and Firm Support Edges for Transfers



Bed Specifications & Dimensions

- 10" Low Height (17" Top of Mattress)
- 32" High Height (39" Top of Mattress)
- Mattress Size: 79"L x 39"/48"W
- Mattress Height: 7" Foam or 8" Air
- Overall Bed Size: 81"L x 44"/52"W
- Weight Capacity: 500 lbs.
- Multi-Height Assist Rail
- 4-Wheel Central Brake Pedal



Available Bed Accessories

- Overbed Reading Light
- Underbed Auto Nightlight
- 8" Length Extension Package
- Additional Set of Rails
- Rail Organizer and Rail Pads
- Overhead Trapeze Helper Bar
- Extra Large Overbed Table
- Sheet Sets & Adjustable Pillows
- Organic Cotton Duvet Set