

Introducing *Supernal* 5



Standard Functions

- Head Adjustable
- Foot Adjustable
- Hi-Low Adjustable
- Trendelenburg
- Reverse Trendelenburg



Quick Facts


- Less than 12 inches from deck to floor
- Height range - 12 inches
- 5 Function Model
- 300 pound weight capacity
- Robust, modular control box design
- Compatible with most Hoyer Lift Systems

Available Sizes

Twin 38" x 80"
Full 53" x 80"
Queen 60" x 80"
Also available in Dual King

Options

- Cherry style finish headboard
- Half Rails
- Bamboo rail covers



Call or email for your free custom evaluation today!
877-445-6233 or sales@transfermaster.com

*Patent Pending