



WHEELCHAIR SEATING SIZE & WEIGHT CHART

7271 Park Circle Drive
Hanover, MD 21076

21055 Ashburn Crossing Dr. Unit 120
Ashburn, VA 20147

www.FreedomMobilitySolutions.com

	Up Thru 5'2" Hemi Height				1. Patient weight and height, look up seat width							
	< 5'0"	5'0"	5'2"	5'4"	5'6"	5'8"	5'10"	6'0"	6'2"	6'4"	6'6"	> 6'6"
<100lbs	14"/16"	14"/16"	14"/16"	16"	16"	16"	16"	16"	16"	16"	16"	16"
100lbs	14"/16"	14"/16"	16"	16"	16"/18"	16"	16"	16"	16"	16"	16"	16"
120lbs	18"	18"	18"	18"	16"/18"	16"/18"	16"	16"	16"	16"	16"	16"
140lbs	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"	18"	16"	16"
160lbs	20"	20"	18"/20"	18"/20"	18"	18"	18"	18"	18"	18"	18"	18"
180lbs	20"	20"	20"	20"	20"	18"/20"	18"/20"	18"	18"	18"	18"	18"
200lbs	22"	22"	20"	20"	20"	20"	20"	18"/20"	18"/20"	18"/20"	18"	18"
220lbs	22"	22"	22"	20"	20"	20"	20"	20"	20"	20"	20"	20"
240lbs	24"	24"	22"/24"	22"/24"	20"	20"	22"	22"	22"	22"	22"	22"
260lbs	24"	24"	24"	24"	22"/24"	22"/24"	22"	22"	22"	22"	22"	22"
280lbs	26"	26"	26"	24"	24"	24"	22"/24"	22"	22"	22"	22"	22"
300lbs	26"	26"	26"	26"	26"	24"	24"	24"	24"	22"	22"	22"
320lbs	28"	28"	28"	26"/28"	26"	26"	24"	24"	24"	24"	24"	24"
340lbs	28"	28"/30"	28"/30"	28"/30"	26"/28"	26"/28"	26"	24"	24"	24"	24"	24"
2. Approx. Seat Length				← 16" →				← 18" →		← 20" →		

Patient under 5'2" Hemi (Lowered) height, 5'4" with 2" cushion
Patient in-between size, go down for men and up for women

This chart is approximate, please measure for EXACT fit!