

# UNIVERSAL HALO SAFETY RING

## PRODUCT GUIDE

PART NUMBER	DESCRIPTION
77301	Halo Assembly with 1 Ring
77302	Halo Assembly with 2 Rings

- \* 77101 has been discontinued and replaced by 77301
- \* 77102 has been discontinued and replaced by 77302
- \* 77103 has been discontinued

### IMPROVEMENTS & UPDATES

The new Halo assembly fits on more beds including all of the Joerns Easy Care Line and Invacare CS Line.

New assembly replaces 77113 and 77203 that were unable to articulate with the bed frame.

Both the 77301 and 77302 come with 2 ring mount brackets and able to mount on either side of the bed for single assemblies.


New assembly hardware includes center bar with mounting holes, screws, and 2 Allen wrenches.

Initial adjustments are based on the bed width and can be done prior to installing on frame using 1 Allen wrench.

Final adjustments are done at each side of the bed frame using 1 Allen wrench.

Also includes cone tip screws that vertically lock into the bed frame to help stabilize and secure the Halo assembly.

The new assembly has improved strength and durability.



### UNIVERSAL HALO FIT LIST

HOSPITAL BED	PART NUMBER
Panacea 1000	77301/77302
Panacea 1500	77401/77402
Panacea 6300 Bariatric	77301/77302
Joerns Care 100 35" W	77301/77302
Joerns UCXT 35" W	77301/77302
Joerns EasyCare Line 35" W	77301/77302
Joerns EasyCare Line 39" W	77301/77302
Joerns EasyCare Line 42" W	77301/77302
Invacare CS Line 36" W	77301/77302
Invacare CS Line 39" W	77301/77302
Invacare CS Line 42" W	77301/77302

HOSPITAL BED	PART NUMBER
Joerns EC2000 (discont.)	77301/77302
Invacare SC900	77301/77302
Invacare 820	77301/77302
Invacare CS C59FX600 36"-42" W	77301/77302
Simmons (All)	77301/77302
Medline 1232	Unable to Attach to Halo Products
Zenith PBH093A35MA	Unable to Attach to Halo Products
Invacare CS9 FX600	Unable to Attach to Halo Products
Noa Elite Riser 35" W	Unable to Attach to Halo Products

\*Universal Halo attaches to additional bed models not listed above.