



Shown in Graphite

Titan PR448



Comfort Zone[™]
by GOLDEN


The Titan is medium-sized and offers a stylish, overstuffed seamed back design, plush chaise pad, and numerous standard features providing hours of unmatched comfort. Your positioning options are endless with our 4-Zone Comfort Zone Positioning System, which includes our patented MaxiComfort[®] Positioning, adjustable headrest and adjustable lumbar all operated by our programmable AutoDrive[™] hand control. Other features include a built-in cup holder, wireless smartphone charger, and accessory mount with a table, conveniently located in the open arm design. Also standard is an additional USB charging port for your phone or table located on the side of the chair.



Standard Fabrics:



Acorn



Graphite



Maple

This faux leather is a rich and modern alternative to natural leather that is extremely comfortable while offering the rugged durability of leather. Made from 100% polyester, the fabric is easy to care for ensuring functionality as well as elegance. This long-lasting fabric will complement a variety of decors while standing the test of time under heavy use.

Crossover fabrics are unavailable.

Standard Features:

- Patented MaxiComfort[®] Positioning
- Adjustable Headrest
- Adjustable Lumbar
- Programmable AutoDrive[™] Hand Control
- USB Charging Port
- Wireless Smartphone Charger
- Accessory Mount with Removable Table
- Built-in Cup Holder
- Extra Pocket
- SmartTek[®] Diagnostic System



DeLUNA Series PR448

Model Number	PR448-MED
Size	Medium
Class II Medical Device	YES
Weight Capacity	375 LBS
Positioning	Twilight
Comfort Zone(s)	4
Twilight Technology	NO
ZG+	NO
Chaise Seat	YES
Pocketed Coil Seat	YES
Back Type	Seam
Removable Back Type	Quick Release
Head Rest Cover	NO
Arm Rest Covers	NO
Overall Width	38.5"
Overall Height	46"
Seat Width	20.5"
Floor to Top of Seat	20.5"
Seat Depth	21"
Seat to Top of Back	31"
Seat to Top of Arm	6.5"
Distance Required from Wall Reclined	26"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	YES
Warranty on Electrical Parts	3 YR
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 YR
Warranty on Electrical Labor	1 YR
Hand Control Routing	Inside Arm
Extra Pocket Option	Standard
Heat & Vibrating Massage Option	NO
HeatWave Technology Option	NO
Articulating Headrest Option	Standard
Power Adjustable Lumbar Option	Standard
Foot Rest Extension Option	NO
Straight Lift Option	Standard
USB Charging Port	Left Arm
Battery Backup	Lithium-ion
SmartTek Diagnostic System	Standard
Place of Origin	Imported



The Titan comes standard with the following comfort and convenience features:



Adjustable Headrest and Lumbar




Accessory Mount with Removable Table



Built-in Cup Holder



Wireless Smartphone Charger



USB Charging Port

Designed and engineered in Old Forge, PA. Assembled in China using Golden components.



Comfort Zone™ by GOLDEN

Comfort Zones 1 - 4 provide an ergonomic comfort experience with an emphasis on head and back support. Each Comfort Zone can be adjusted independently by using the exclusive AutoDrive hand control.



LOWER BODY


UPPER BODY

HEADREST

LUMBAR

TWILIGHT TILT

WE PATENTED COMFORT™



Scan with your smartphone camera to learn more about Comfort Zone by Golden.



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Complies with
UL 73 and
CSA C22.2 No. 68



Find us on Facebook
www.facebook.com/goldentechnologies

SHT-PR448-011124