

Full-Electric Bariatric Bed, 42"

HCPCS: E0302/E0304/E1399



FIGURE A

- Quiet, smooth operation
- Heavy-duty frame ensures strength and patient safety
- All steel construction
- Larger sleep surface than conventional bed
- Split-pan design with removable bed ends is easy to set up
- No tools necessary for set up
- Hand control (included) provides for multiple bed positioning for patients
- Low air loss/alternating pressure mattress available
- "T" rails available (Sold separately) (Figure A)
- Emergency manual crank included
- Total Weight Capacity 600 lbs, Maximum Patient Weight 500 lbs.

PRODUCTS

	Description
15300	Bariatric Bed, Full Electric, 42", 1/cs
15300BV-1HR	Bariatric bed w/T rail, 1pkg/cs
15300BV-PKG	Bariatric Bed Pkg, 1pkg/cs

SPECIFICATIONS

	Description
Bed Height	18.25" - 26.25"
Carton Shipping Weight	276 lbs.
Sleeping Surface	80"
Unit Dimensions	42" (W) x 88" (L)
Warranty	(Original Purchaser) Limited Lifetime on Welds, 5 Years on Frame, 1 Year on All Other Parts and Components
Weight Capacity	600 lbs
Bed Min height	18.25"
Bed Max height	26.25"

ACCESSORIES

	Description
15301	Bariatric Mattress, 42in, 1/ea, EA

PARTS

	Description
15307	Head&Foot Board Set, 15300-PKG, 1/cs, CS
15003HC	Head Crank Assembly 15003, 1ea, EA
15300P1001	Head Motor Ext. Cord, Bariatric, 1/ea, EA
15005SPN	Bed HiLo Shaft Pin, New, 1ea, EA
15044CNL	Caster, NoLock for 15544, 15547, 1/ea, EA

PARTS	Description
15044CL	Caster,w/Lockfor15544,15547,1/ea, EA
15005PF1028	Gearbox,Hi/LowShaft,FootBrd,1/ea, EA
15005P-1040	Hooks to Attch Sleep Surface,6/bx, BX
15005-SPRING	HelicalSpring/SleepSurf,Top&Bottom,1/ea, EA
15005P1047	Foot Crank Handle,1/ea, EA
15005P1045	CrankHandle,Man.Bed,1/ea, EA
15005PS1051	Panel Frame Screw,Bed,1/ea, EA
15005P1046	FootCrank Gear w/Handle,1/ea, EA
15005P1028	GearBox,Hi/LowShaft,HeadBrd,1/ea, EA
15005S	Hi/Low Shaft forSemi/Full,15300,1/ea, EA
15015F	Foot Board,SemiElec.Bed,1/cs, CS
15015H	Head Board,Full&SemiBed,1/cs, CS
15008	Head Spring,Man,Semi,Full,1/cs, CS
15005P1051	Panel Frame for Headboard,1/ea, EA
15005P1151	Panel Frame for Footboard,1/ea, EA