

Tuffy Wide Recliner™ 497

Extra Wide Recliner

SEAT WIDTH AVAILABLE: 22", 24"
UPHOLSTERY AVAILABLE: DARK BLUE



The Tuffy Recliner™ offers both maximum durability and comfort for larger individuals who seek long term rehabilitation. With its heavy duty frame and precision engineered hydraulic assist reclining mechanism, this chair is ready for frequent use without any problem.

Tuffcare®
HIGHER STANDARD FOR HEALTH CARE™

Tuffy Wide Recliner™ 497

Extra Wide Recliner

SEAT WIDTH AVAILABLE: 22", 24"
UPHOLSTERY AVAILABLE: DARK BLUE

STANDARD FEATURES

- * Folding triple chrome plated frame with detachable desk arms
- * Hydraulic assist cylinders recline back from 90% to 180%
- * Padded double embossed, flame retardant vinyl upholstery, padded arm pads & calf pads
- * 8" x 1 3/4" cast wheels with pneumatic airless tire
- * 24" x 1 3/8" rear molded wheels with pneumatic airless tire
- * Adjustable toggle wheel locks
- * Cam action swing away detachable elevating legrests
- * Stainless side panels & aluminum foot plates
- * Weight limit: 22" 370 lbs, 24" 400 lbs


OPTIONAL ACCESSORIES

- * Anti-theft bar
- * Wheelchair tray
- * Anti-folding device
- * Seat belt
- * Projection handrim
- * Anti-theft lock for front riggings
- * Pneumatic tires
- * I.V. rod attachment
- * Anti-tipping device
- * Headrest extension
- * Wheel lock extension tubes
- * Leg straps
- * Oxygen tank holder


SPECIFICATIONS

MODEL	SEAT			BACK	FRONT RIGGINGS	OVERALL					WEIGHT (w/o Front Riggings)
	22"	20"	19"			30"	42"	50 1/2"	37"	13"	
4979E	22"	20"	19"	33"	16"t 21"	30"	42"	50 1/2"	37"	13"	56 lbs.
497XE	24"	20"	19"	33"	16"- 21"	32"	42"	50 1/2"	37"	13"	58 lbs.

E ELEVATING LEGRESTS N SR. EXTRA WIDE SEAT 24" WIDTH



Reinforcement weld on all key locations and double cross bar construction makes this the most reliable heavy-duty recliner in the market.



Standard with front and back airless pneumatic tires for wider track and durability.



With heavy-duty hydraulic assist cylinder, lifting heavy patient from a reclining position never been easier.