

Fall 2017
nova

GO!
MOBILITY





NOVA Transport Chairs and Wheelchairs **bring people together** with the best features and quality in the mobility world.

FRAME STRENGTH

NOVA Transport Chairs and Wheelchairs are made with the highest quality materials and workmanship. They are backed by our Limited Lifetime Warranty.

BEST BRAKES

Whether it is the locking hand brakes on some models or wheel locks on others - NOVA has the best brakes around!

UPHOLSTERY

You can see and definitely feel the difference with NOVA's seating upholstery. The material is attractive and comfortable, the craftsmanship is excellent, and the quality will last through life's many adventures.

FOOTRESTS

All NOVA Transport Chairs and Wheelchairs come with swing away, removable footrests. They are easy to adjust and durable for resting the unique shapes and sizes of our feet.

THE EXTRAS

With a NOVA, you get the Extras as standard that come at an additional cost with other brands. Our 6000 series Reclining and 7000 series Lightweight Wheelchairs come standard with Anti-Tippers and all Transport Chairs come standard with seat belts.

Tips for choosing your NOVA Transport Chair or Wheelchair

What's the main difference between a Wheelchair and a Transport Chair?

Wheelchairs can be self-propelled by the user and also pushed by the companion. The large 24" rear wheels with rims allow the user to self-propel. This feature does add more weight and bulk to the wheelchair. Transport chairs have smaller rear wheels (8" & 12") so the user cannot self-propel and can only be pushed by the companion. Thus, the transport chair is both lighter and more compact than a wheelchair.

The Brakes - Who's in Control? (Option for Transport Chairs)

Some NOVA Transport Chairs come with both rear wheel locks and locking hand brakes at the companion handles. This allows the companion to control both the speed of the ride and ability to park the brake. This is a great option (& promotes safety) for those bumpy or hilly terrains. Make sure to have your seat belt on!

Weight and Width Matter!

Be sure to understand and match the weight capacity and the seat width of the transport chair or wheelchair to the weight and general size of the user. You should also carefully consider the actual weight and size of the wheelchair or transport chair, especially if there are lots of events and adventures in the plans.

RIDE IN STYLE

Choose the color and options that compliment your own personal tastes and lifestyle. Add a Leopard Mobility Clutch for your cell phone, a Flashlight Attachment for the evening, a Butterfly Mobility Bag for your important belongings, and a Cup Holder for your favorite beverages!

 *A great value!*



COMES STANDARD

Carbon steel hammertone frame
Fixed, full-length arms
Quick release fold-down back
Locking rear wheels
Non-skid foot plates & heel loops
Foldable

OPTIONS

Travel Bag
Cup Holder
Flashlight attachment
Designer Mobility Bags
Color: Hammertone

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19"

WIDTH BETWEEN ARMS: 16.5"

ARM HEIGHT FROM SEAT: 11.25"

BACK HEIGHT FROM SEAT: 18.25"

WHEEL SIZE: 8"

OVERALL DIM: 39.5"h x 23.5"w x 27"d

DIM WHEN FOLDED: 31.5"h x 10"w x 26.5"d

WEIGHT: 23 lb

SEAT DIM: 19"w x 15.75"d



 *A great value!*



COMES STANDARD

Carbon steel frame
Fixed, full-length arms
Quick release fold-down back
Locking rear wheels
Non-skid foot plates & heel loops
Foldable

OPTIONS

Travel Bag
Cup Holder
Flashlight attachment
Designer Mobility Bags
Colors: Black, Blue, Champagne, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19.75"

WIDTH BETWEEN ARMS: 16.5"

ARM HEIGHT FROM SEAT: 11"

BACK HEIGHT FROM SEAT: 17.75"

WHEEL SIZE: 8"

OVERALL DIM: 39.75"h x 23.5"w x 27.5"d

DIM WHEN FOLDED: 31.5"h x 9.75"w x 26.5"d

WEIGHT: 22.75 lb

SEAT DIM: 19"w x 15.75"d



✓ *Lightweight!*



COMES STANDARD

- Lightweight aluminum frame
- Fixed, full-length arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- 327 Colors: Blue, Red
- 329 Colors: Black, Blue, Red

WEIGHT CAPACITY: 300 lb

WEIGHT: 19.5 lb

SEAT HEIGHT: 19.5"

SEAT DIM: (327) 16.75" w x 15.75" d
(329) 19" w x 15.75" d

WIDTH BETWEEN ARMS: (327) 14.5" & (329) 16.5"

ARM HEIGHT FROM SEAT: 11"

BACK HEIGHT FROM SEAT: 18.25"

WHEEL SIZE: 8"

OVERALL DIM: (327) 39.5" h x 20.75" w x 28" d
(329) 39.5" h x 23.5" w x 28" d

DIM WHEN FOLDED: 31.25" h x 9.5" w x 26.75" d



✓ *Lightweight!*



COMES STANDARD

- Lightweight aluminum frame
- Fixed, full-length arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Tan with Plaid Upholstery

WEIGHT CAPACITY: 300 lb

WEIGHT: 19.5 lb

SEAT HEIGHT: 19.5"

SEAT DIM: 19" w x 15.75" d

WIDTH BETWEEN ARMS: 16.5"

ARM HEIGHT FROM SEAT: 11"

BACK HEIGHT FROM SEAT: 18.25"

WHEEL SIZE: 8"

OVERALL DIM: 39.5" h x 23.5" w x 28" d

DIM WHEN FOLDED: 31.25" h x 9.5" w x 26.75" d



 Desk Arms are great for dining!



COMES STANDARD

- Lightweight aluminum frame
- Fixed desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Blue

WEIGHT CAPACITY: 250 lb

SEAT HEIGHT: 19"

WIDTH BETWEEN ARMS: 17.5"

ARM HEIGHT FROM SEAT: 8.75"

BACK HEIGHT FROM SEAT: 15.25"

WHEEL SIZE: 8"

OVERALL DIM: 37"h x 22.5"w x 29.5"d

DIM WHEN FOLDED: 29.5"h x 9.25"w x 23.25"d

WEIGHT: 20 lb

SEAT DIM: 17.75"w x 16"d



 Desk Arms are great for dining!



COMES STANDARD

- Lightweight aluminum frame
- Fixed desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Black, Blue, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19"

WIDTH BETWEEN ARMS: 17.5"

ARM HEIGHT FROM SEAT: 9"

BACK HEIGHT FROM SEAT: 15.25"

WHEEL SIZE: 8"

OVERALL DIM: 37"h x 22.5"w x 30"d

DIM WHEN FOLDED: 29.5"h x 9.75"w x 23.75"d

WEIGHT: 18.75 lb

SEAT DIM: 18"w x 16"d



✓ *Hand Brakes and 12" wheels.
Go be adventurous!*



COMES STANDARD

- Lightweight aluminum frame
- Hand brakes for safety/control
- Fixed, full-length arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19.5"

WIDTH BETWEEN ARMS: 17.25"

ARM HEIGHT FROM SEAT: 12.5"

BACK HEIGHT FROM SEAT: 19.25"

WHEEL SIZE: 12" rear wheels, 8" front wheels

OVERALL DIM: 40.5"h x 24"w x 30.5"d

DIM WHEN FOLDED: 32.5"h x 10"w x 29.25"d

WEIGHT: 24.75 lb

SEAT DIM: 19.75"w x 15.75"d



✓ *The ultimate ... Brakes, Big Wheels, and
Flip-Up Arms!*



COMES STANDARD

- Lightweight aluminum frame
- Hand brakes for safety/control
- Flip-up desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19.75"

WIDTH BETWEEN ARMS: 19.25"

ARM HEIGHT FROM SEAT: 10"

BACK HEIGHT FROM SEAT: 17.75"

WHEEL SIZE: 12" rear wheels, 8" front wheels

OVERALL DIM: 38.25"h x 25.5"w x 31"d

DIM WHEN FOLDED: 30.5"h x 12.5"w x 29.5"d

WEIGHT: 27.5 lb

SEAT DIM: 18.5"w x 15.5"d



 Move over easily with Flip-Up Arms!



COMES STANDARD

- Lightweight aluminum frame
- Flip-up desk arms
- Quick release fold-down back
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 300 lb

SEAT HEIGHT: 19"

WIDTH BETWEEN ARMS: 19.25"

ARM HEIGHT FROM SEAT: 10"

BACK HEIGHT FROM SEAT: 18"

WHEEL SIZE: 8"

OVERALL DIM: 38.5"h x 25"w x 28.5"

DIM WHEN FOLDED: 30.5"h x 12.5"w x 26.75"d

WEIGHT: 23 lb

SEAT DIM: 18.5"w x 15.5"d



 Extra-Wide and Extra-Strong!



COMES STANDARD

- Lightweight aluminum frame
- Hand brakes for safety/control
- Flip-up desk arms
- Locking rear wheels
- Non-skid foot plates & heel loops
- Foldable

OPTIONS

- Travel Bag
- Cup Holder
- Flashlight attachment
- Designer Mobility Bags
- Colors: Blue, Red

WEIGHT CAPACITY: 400 lb

SEAT HEIGHT: 20"

WIDTH BETWEEN ARMS: 22"

ARM HEIGHT FROM SEAT: 9.75"

BACK HEIGHT FROM SEAT: 17.75"

WHEEL SIZE: 12" rear wheels, 8" front wheels

OVERALL DIM: 38"h x 28.5"w x 30"d

DIM WHEN FOLDED: 38"h x 11"w x 30.5"d



WEIGHT: 27.75 lb

SEAT DIM: 21.5"w x 15.5"d

K1



16"/18" WIDTHS

5060S/5080S - Fixed full arms with easy adjust footrests
 5060SE/5080SE - Fixed full arms with elevating leg rests

FEATURES

Hammertone finish
 Black vinyl upholstery
 Dual axle positions for quick seat height conversion

WEIGHT CAPACITY: 300 lb

WEIGHT: (5060S) 34 lb / (5080S) 36 lb

SEAT HEIGHT: 20"

SEAT DIM: (5060S) 16"w x 16"d / (5080S) 18"w x 16"d

WIDTH BETWEEN ARMS: (5060S) 13.75" / (5080S) 15.75"

ARM HEIGHT FROM SEAT: (5060S) 9.5" / (5080S) 10"

BACK HEIGHT FROM SEAT: (5060S) 15" / (5080S) 14.75"

OVERALL DIM: (5060S) 36"h x 23.25"w x 30.25"d
 (5080S) 36"h x 24.5"w x 31.25"d

FOLDED WIDTH: (5060S) 10.75" / (5080S) 11.25"



K1



16"/18" WIDTHS

5165S/5185S - Detachable desk arms with easy adjust footrests
 5165SE/5185SE - Detachable desk arms with elevating leg rests
 5166S/5186S - Detachable full arms with easy adjust footrests
 5166SE/5186SE - Detachable full arms with elevating leg rests

FEATURES

Hammertone finish
 Black vinyl upholstery
 Dual axle positions for quick seat height conversion

WEIGHT CAPACITY: 300 lb

WEIGHT: (5165S) 38 lb / (5185S) 39 lb

SEAT HEIGHT: 20"

SEAT DIM: (5165S) 16"w x 16"d / (5185S) 18"w x 16"d

WIDTH BETWEEN ARMS: (5165S) 16.25" / (5185S) 18.25"

ARM HEIGHT FROM SEAT: 9"

BACK HEIGHT FROM SEAT: 15.25"

OVERALL DIM: (5165S) 35.25"h x 24"w x 31.5"d
 (5185S) 35.25"h x 26"w x 31.5"d

FOLDED WIDTH: (5165S) 12" / (5185S) 12.25"



K3



16"/18" WIDTHS

- 5160S/5180S - Detachable desk arms with easy adjust footrests
- 5160SE/5180SE - Detachable desk arms with elevating leg rests
- 5161S/5181S - Detachable full arms with easy adjust footrests
- 5161SE/5181SE - Detachable full arms with elevating leg rests

FEATURES

- Hammertone finish
- Black vinyl upholstery
- Dual axle positions for quick seat height conversion

WEIGHT CAPACITY: 300 lb

WEIGHT: 33 lb

SEAT HEIGHT: 20"

SEAT DIM: (5160S) 16"w x 16"d / (5180S) 18"w x 16"d

WIDTH BETWEEN ARMS: (5160S) 16.25" / (5180S) 18.25"

ARM HEIGHT FROM SEAT: 9"

BACK HEIGHT FROM SEAT: 15.25"

OVERALL DIM: (5160S) 35.5"h x 24.25"w x 31.5"d
(5180S) 35.5"h x 26.25"w x 31.5"d

FOLDED WIDTH: 12.75"



K7



20" WIDTH

- 5200S - Detachable desk arms with easy adjust footrests
- 5200SE - Detachable desk arms with elevating leg rests
- 5201S - Detachable full arms with easy adjust footrests
- 5201SE - Detachable full arms with elevating leg rests

FEATURES

- Hammertone finish
- Black vinyl upholstery
- Dual axle positions for quick seat height conversion

WEIGHT CAPACITY: 300 lb

WEIGHT: 33.5 lb

SEAT HEIGHT: 20"

SEAT DIM: 20"w x 16"d

WIDTH BETWEEN ARMS: 20.25"

ARM HEIGHT FROM SEAT: 9"

BACK HEIGHT FROM SEAT: 15.25"

OVERALL DIM: 35.5"h x 28.25"w x 31.5"d

FOLDED WIDTH: 13"



K7



22"/24" WIDTHS

- 5220S/5240S - Detachable desk arms with easy adjust footrests
- 5220SE/5240SE - Detachable desk arms with elevating leg rests
- 5221S/5241S - Detachable full arms with easy adjust footrests
- 5221SE/5241SE - Detachable full arms with elevating leg rests

FEATURES

- Hammertone finish
- Black vinyl upholstery
- Dual axle positions for quick seat height conversion
- Dual cross bar supports for added strength

WEIGHT CAPACITY: (5220S) 350 lb / (5240S) 425 lb

WEIGHT: (5220S) 53 lb / (5240S) 54 lb

SEAT HEIGHT: 20"

SEAT DIM: (5220S) 22"w x 18"d / (5240S) 24"w x 18"d

WIDTH BETWEEN ARMS: (5220S) 22.75" / (5240S) 24.75"

ARM HEIGHT FROM SEAT: (5220S) 9" / (5240S) 8.5"

BACK HEIGHT FROM SEAT: 15.25"

OVERALL DIM: (5220S) 36"h x 30.5"w x 33.75"d
(5240S) 36"h x 32.5"w x 33.75"d

FOLDED WIDTH: (5220S) 12.75" / (5240S) 13"



K4



16"/18" WIDTHS

- 7160L/7180L - Flip back desk arms with easy adjust footrests
- 7160LE/7180LE - Flip back desk arms with elevating leg rests
- 7161L/7181L - Flip back full arms with easy adjust footrests
- 7161LE/7181LE - Flip back full arms with elevating leg rests

FEATURES

- Hammertone finish
- Black nylon upholstery
- Dual axle positions for quick seat height conversion

WEIGHT CAPACITY: 300 lb

WEIGHT: 32 lb (7160L) / 33 lb (7180L)

SEAT HEIGHT: 20"

SEAT DIM: (7160L) 16"w x 16"d / (7180L) 18"w x 16"d

WIDTH BETWEEN ARMS: (7160L) 16.25" / (7180L) 17.75"

ARM HEIGHT FROM SEAT: 10.5"

BACK HEIGHT FROM SEAT: 16.25"

OVERALL DIM: (7160L) 36.25"h x 23.5"w x 31.5"d
(7180L) 36.25"h x 25.5"w x 31.5"d

FOLDED WIDTH: 11.75"



K4



20" WIDTH

7200L - Flip back desk arms with easy adjust footrests

7200LE - Flip back desk arms w/ elevating leg rests

7201L - Flip back full arms with easy adjust footrests

7201LE - Flip back full arms w/ elevating leg rests

FEATURES

Hammertone finish

Black nylon upholstery

Dual axle positions for quick seat height conversion

Dual cross bar supports for added strength

WEIGHT CAPACITY: 300 lb

WEIGHT: 33 lb

SEAT HEIGHT: 20"

SEAT DIM: 20" w x 16" d

WIDTH BETWEEN ARMS: 20.25"

ARM HEIGHT FROM SEAT: 10.5"

BACK HEIGHT FROM SEAT: 16.25"

OVERALL DIM: 36.25" h x 27.5" w x 31.5" d

FOLDED WIDTH: 11.75"



MOBILITY BAGS



4009BK 4008GA 4008LP 4008BB 4008EG

4008GA (Enchanted Garden)

4008LP (Leopard Print)

4008BB (Boho Blossoms)

4008EG (English Garden)

4009BK (Classic Black)

HANGING MOBILITY POUCH

4001WP (Black)

P4200SB (Aloha)

- Convenient for carrying various personal belongings
- Easily attaches to all Wheelchairs, Transport Chairs, and Walkers



P4200SB

ANTI-TIPPERS

- Helps prevent wheelchairs from tipping backward



Item #	Fits/Serial Numbers/Color
W-5105	5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome
W-5106	5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome
W-5005	5200S,5180S,7160S,7200L,7180L,7160L/Y,J,H/Black
J-AT101	5060S,5080S,5165S,5185S,5220S,5240S/Y,J,H/Chrome

ELEVATING LEGRESTS

W-7001

- For 5000/6000/7000 series



The NOVA Promise




Over twenty years ago, we made the commitment to put our customers first. We honor our warranties and we care. Simply put - We will do whatever it takes.

NOVA is a leading innovator and manufacturer of Mobility, Bathroom, and Independent Living products including canes, rolling walkers, transport chairs, and mobility accessories. Since 1993, NOVA has grown organically with a strong foundation and values of service, quality and integrity. NOVA's mission is to provide superior products with great function and fashion for people to live a healthy, independent, and beautiful lifestyle. Our commitment is represented by our extraordinary staff, excellent warranty, education, innovation, and charitable contributions globally.

NOVA's commitment to customers begins by manufacturing each product with the highest quality standards. Our commitment continues for the life of the product with caring and professional customer support. **All NOVA Transport Chairs and Wheelchairs are backed by a Limited Lifetime Warranty.**

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

CARE • QUALITY • INNOVATION




THE NOVA PROMISE

NOVA Transport Chairs and Wheelchairs are **Built for Life** with the best features and quality in the mobility world.



Learn More Online

To learn more about: choosing the perfect rolling walker for you, assembling your walker, product features, safety tips, accessories and so much more, watch the latest videos on:

-  (NOVA Medical Products)
- NOVA Website: novajoy.com



Connect with NOVA!

novajoy.com
facebook: [novamedicalprod](https://www.facebook.com/novamedicalprod)
twitter: [novamedicalprod](https://twitter.com/novamedicalprod)
youtube: [NOVA Medical Products](https://www.youtube.com/NOVAMedicalProducts)





**LIVE LIFE TO THE
FULLEST WITH YOUR NOVA:
The very best in Transport Chairs
and Wheelchairs.**

We are passionately committed to connecting human relationships and fulfilling human experiences by providing products and services that help fully realize the potential and possibilities of people with physical challenges. NOVA is invested in the Human Ability through our Products, Services and Culture.

WARNING

- Lock wheel locks when sitting stationary
- Do not use on stairs, escalators, or moving walkways
- Consult your user manual before operating



Authorized NOVA Dealer