

TRILLIUM Series



THE TRILLIUM SERIES IS THE MOST COMFORTABLE WALKER ON THE MARKET TODAY.

With the new padded seat cushion and wide back strap pad, users can sit in comfort and walk in style. The Trillium comes with an all-new braking system, allowing for effortless handle height adjustment and easy maintenance.

FEATURES

SEAT 1.5" thick padded seat cushion for extra comfort

BRAKE SYSTEM Concealed cable / Easy to adjust and maintain

COLLAPSIBLE BASKET Fully covered for increased privacy. Foldable & removable

WHEELS 8" lightweight polyurethane tires. Offset rear forks for wider gait

WIDE BACK STRAP PAD Extra wide for added comfort

HAND GRIPS Removable design for easy maintenance

COLORS Blue, Carbon Black, Pink, Shiraz Red



SPECIFICATIONS

MINI

14 lbs*
18" Seat Width
18" Seat Height
26" Width x 25.5" Length
29" - 32" Handle Height
300 lbs Weight Capacity

REGULAR

14 lbs*
18" Seat Width
21" Seat Height
26" Width x 25.5" Length
31.5" - 34.5" Handle Height (optional tall handles)
34.5" - 37.5" in Carbon Black)
300 lbs Weight Capacity

TALL

14 lbs*
18" Seat Width
24" Seat Height
26" Width x 25.5" Length
34.5" - 37.5" Handle Height (optional tall handles)
37.5" - 40" in Carbon Black)
350 lbs Weight Capacity

OPTIONAL ACCESSORIES

- Cup Holder
- Cane Holder
- Phone Holder
- Slow Down Brakes
- Eyeglass Holder
- Tray
- Curb Climber

*Net Weight of walker without accessories, backstrap and folding basket



Evolution
TECHNOLOGIES

Toll Free 1 800 556 2558

info@evolution.ca

www.evolution.ca