




Standard Colors



White



Teal



Black

**Remote
Control App
for Enhanced
Convenience**

**Folds Easily
Drives Remotely
Weighs Only 56 lbs.**

**Ally
GP303**



Visit App Store
for Remote
Control



Experience freedom with the new Golden Ally™, our new power chair designed to provide the strength and solid feel of a traditional chair yet remaining light weight and portable enough to be the perfect travel partner. Enjoy the maneuverability indoors with its tight turning radius or take on an outdoor adventure and appreciate its great ride quality and capability, compliments of its large, pneumatic rear drive wheels. Its comfortable seating and foot positions will be a welcome change from the compromises of the typical travel chair. With a long-lasting battery, user-friendly interface, and Smart Control Mobile App, navigate any environment with ease and confidence. Move effortlessly with the Ally – your perfect travel partner!

Standard Features:

- Simple Folding Mechanism
- Smart Control Mobile App
- Tight Turning Radius
- Memory Foam Seat Cushion
- Ergonomically Designed Seating & Foot Positions
- Airline Friendly Lithium-Ion Battery
- System Monitoring Voice Alerts
- A&C USB Ports
- Front & Rear Lighting



Model Number	GP303
Class II Medical Device	Yes
Medicare Code	Not Coded
Weight Capacity	330 lbs.
Drive Wheels	Rear
Maximum Speed	5 mph
Operating Range ¹	9 miles
Ground Clearance	2.5"
Turning Radius	28.5"
Type Batteries	12 AH Lithium-Ion
Colors	Green, White, Black
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp Remote
Dimensions:	
Length	41"
Width	24.5"
Width Between Arms	18.5"
Height (ground to top of back of seat)	36"
Ground to Top of Seat	21.5"
Weight of Unit: (assembled)	
Heaviest Piece	49.5 lbs.
Battery	5 lbs.
Tires:	
Drive Wheels	10" Pneumatic
Casters	8" Solid
Anti-Tip Wheels	2"
Standard Seat:	
Back Height - no headrest	18"
Width x Depth	17" x 16"
Color	Black Fabric
Footrest:	
Height adjustable	Yes
Increments	1"
Warranty:	
Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months
Batteries	12 Months

¹ Battery range tested at 200 lbs. and will vary due to rider weight, drive surface, terrain and battery type.

² Electronics warranty excludes batteries.



Zippered Waterproof Backpack



Highly Visible Rear Lighting



Controller is pre-wired allowing right or left side placement

Optional Accessories:



Extendable Headrest



Cup/Phone Holder



Heated Seat Control



Wireless Remote



Phone: 800-624-6374

Fax: 800-628-5165

www.goldentech.com



www.facebook.com/goldentechnologies

SHT-GP303-081524