## WORLD'S LIGHTEST SCOOTER!

## Feather SCOOTER 27X

At just 27 lbs.\*, the Feather<sup>™</sup> Scooter 27X is the lightest mobility scooter on the planet.

Easily take the **NEW Feather™ Scooter 27X** with you everywhere! Weighing **only 27 pounds\*** and **folding** to a compact 16" x 24" - this scooter is the **perfect travel companion**. Its compact size makes it easy to store and is **cruise line compliant** to fit under a cabin bed.

## FEATURES:

- only 27 lbs\*
- 6 mile drive range
- one-fold frame
- airline compliant battery
- cruise line ready

## **SPECIFICATIONS**

| Weight Capacity                | 250 lbs. | Ground Clearance                            | 2"                 |
|--------------------------------|----------|---|--------------------|
| Product Weight w/out batteries | 27 lbs.  | Front Wheel Size                            | 6"                 |
| Product Weight w/batteries     | 30 lbs.  | Rear Wheel Size                             | 6"                 |
| Seat Width                     | 16"      | Turning Radius                              | 51"                |
| Seat Depth                     | 11"      | Max Climbing Angle                          | 6 degrees          |
| Seat-to-Floor Height           | 19"      | Max Speed                                   | 4 mph              |
| Overall Width                  | 16"      | Driving Range                               | 6 Miles            |
| Overall Length                 | 24"      | Battery Type                                | 24V 6.6 AH         |
| Overall Height                 | 32"      | Warranty on Frame                           | Three Year Limited |
| Folded Depth                   | 16"      | Warranty on Wearable<br>Parts & Electronics | One Year Limited   |
| Folded Height                  | 10"      | Warranty on Batteries                       | Six Months Limited |

weight without battery. ^removable lithium battery is airline compliant. Check with individual carrier for guideline

Lightweight

Foldable





Locking Brakes Arrives Assembled





www.FeatherMobility.com