



ProBasics K4 Wheelchair

HIGH STRENGTH LIGHTWEIGHT WHEELCHAIR

The **ProBasics K4 Wheelchair** is a high-strength lightweight extremely versatile wheelchair. With an included seat and upholstery extension combined with all of its other unique features that come standard, the ProBasics K4 easily adapts to any patient's individual needs. It offers superior value and functionality for those looking for a high-performance rental product at a moderate price.



Flip-back armrests are height adjustable



2" upholstery seat extension



Adjustable anti-tippers with push-pin flip-up wheels



Quick release rear axle

FEATURES

- Quick-release axles allow for tool-free removal of rear wheels for fast hemi-height adjustment
- Comes included with a seat rail and upholstery extension that easily adjusts seat depth from 16" to 18" without any additional upholstery
- Flip-back padded arm rests are height adjustable from 10" to 14" to accommodate a seat cushion
- Adjustable height back and adjustable angle back from 5° - 20°
- Comfortable nylon upholstery resists mildew and bacteria
- Chart pocket on backrest for additional convenience
- Includes tool-free adjustable front riggings with composite angle adjustable footplates for more relaxed foot positioning
- Elevating padded leg rests are equipped with a pneumatic piston for infinite elevating positions
- Adjustable anti-tippers with push-pin flip-up wheels
- Easy push-to-lock wheel locks
- HCPCS* Code: K0004

K0195 or E0990 (elevating legrests only)

E0971 (anti-tippers)

K0040 (angle adjustment footplates)

E2201 (non-standard seat frame - 20" chair only)

SPECIFICATIONS

Seat Width:	16" 18" 20"
Seat Depth:	16" - 18"
Seat-to-Floor Height ¹ :	Adult 19.5" Hemi 17.5"
Back Height:	16" - 18"
Back Style:	Height Adjustable
Overall Width w/Arms:	24.5" 26.5" 28.5"
Overall Depth:	32" (without riggings)
Weight Capacity:	300 lb
Product Weight ² :	32 - 36 lb
Carton Dimensions:	33" x 13" x 37"
Shipping Weight:	35 - 39 lb
Warranty:	Limited lifetime on frame One year on all parts

Item #	Seat	Arm Style	Leg Style	UPC	UOM
WC41616DS	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070756	each
WC41816DS	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070763	each
WC42016DS	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Swing-Away Footrests	815067070770	each
WC41616DE	16" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070787	each
WC41816DE	18" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070794	each
WC42016DE	20" x 16"	Height Adjustable, Flip Back, Desk-Length	Elevating Legrests	815067070800	each

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.

¹NOTE: These heights can vary ¼-inch due to tire wear and manufacturing tolerances.

²NOTE: Weights based on 18" x 16" wheelchair with desk-length arms and no front riggings. Weights will vary depending on how the wheelchair is equipped.