



A

B

C

## Poly-Fly High Strength, Lightweight Wheelchair/Flyweight Transport Chair Combo

Item #	Seat Width	Color	UOM
TR16	16"	Silver Frame and Black Upholstery	1/cs
TR18	18"	Silver Frame and Black Upholstery	1/cs
TR20	20"	Silver Frame and Black Upholstery	1/cs

- Can be used as standard self-propelled wheelchair or a transport chair all in one.
- Quick-release 24" wheels can be removed to transition self-propelled chair to a transport chair. (Figure A - B)

### specifications

Overall Width with 24" Wheels (TR16)	23.5" (open)   12"(closed)
Overall Width with 24" Wheels (TR18)	24.5" (open)   12"(closed)
Overall Width with 24" Wheels (TR20)	25.5" (open)   12"(closed)
Overall Height w/ 24" Wheels	37"
Seat-To-Floor Height	19.5"
Back Height	16"
Seat Depth	16"
Overall Length w/ Riggings	42"
Weight w/ 24" Wheels (TR16   TR18   TR20)	33 lbs.   34 lbs   36 lbs.
Carton Shipping Weight w/ 24" Wheels (TR16   TR18   TR20)	38 lbs.   39 lbs.   41 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime

- All-aluminum frame.
- Comes with two sets of wheel locks: One for use with wheelchair, one for use with transport chair.
- Back folds down with deluxe back release. (Figure C)
- Nylon upholstery.
- Comes standard with swing-away footrests with composite foot plates and heel loops.
- 8" casters.
- Comes with anti-tippers.
- STDELR-TF - Swing-away Elevated Leg rest, Tool-free, (Optional front riggings).
- Weight Capacity: 250 lbs.