



# Venture- Euro Style Rollator

**300 LBS**  
WEIGHT CAPACITY



ITEM #	COLOR
810AS	AMBER SUNRISE
810BL	BLUE
810GG	GRAPHITE GREY
810RD	RED

## FEATURES

- Sleek euro style design
- Lightweight aluminum
- Convenient storage bag included
- Cane holder included
- 10" front wheels for indoor and outdoor use

## SPECIFICATIONS

Weight:	15 lbs
Handle Height:	33" - 38"
Seat to Floor Height:	21.5"
Seat Dimensions:	18.75" W x 10.5" D
Unit Dimensions:	24.5" x 27.5"
Weight Capacity:	300 lbs
Available Color(s):	Amber Sunrise, Blue, Graphite Grey, & Red

## WARRANTY

Limited Lifetime warranty against manufacturing defects.

This warranty does not cover product failure due to misuse, negligence, improper storage or handling, improper operation or unauthorized modifications.

For warranty service, it is recommended but not required that the product be returned to the Lifestyle Mobility Aids dealer through whom it was originally purchased. However, all Lifestyle Mobility Aids Dealers are qualified to assist you in obtaining warranty service. If the product is to be returned to Lifestyle Mobility Aids, prior authorization will be required. In the event of a defect covered by this warranty, Lifestyle Mobility Aids will determine if the product will be repaired or replaced.

This Lifestyle Mobility Aids product is identified by a serial number at the time of manufacture. Removing or defacing the serial number label may void this warranty. If you have questions about your Lifestyle Mobility Aids product or this warranty, please contact us at [contactus@lifestylemobilityaids.com](mailto:contactus@lifestylemobilityaids.com).



LIFESTYLE MOBILITY AIDS, INC  
 Email : [contactus@lifestylemobilityaids.com](mailto:contactus@lifestylemobilityaids.com)  
 Website : [www.lifestylemobilityaids.com](http://www.lifestylemobilityaids.com)  
 Phone : 877-843-6464

## ADDITIONAL PRODUCTS TO CONSIDER:



TRANSFER  
BENCH



RAISED HINGED  
TOILET SEAT



TRANSPORTER-  
ALUMINUM WHEELCHAIR/  
TRANSPORT CHAIR