



RIVAL®  
shown in *Evolution Orange*  
with TRU-Balance® 3  
Power Positioning

# **RIVAL**®

The Rival® rear-wheel drive power chair features innovative design that provides exceptional outdoor performance and handling while delivering outstanding tight space maneuverability.

The Rival accepts a complete range of rehab seating and electronics options.

**QUANTUM**®

# Specifications

<b>DRIVE WHEELS</b>	<b>Standard:</b> 14" Solid <b>Optional:</b> 14" Pneumatic; 14" Knobby
<b>CASTER WHEELS</b>	<b>Front:</b> 9" Pneumatic
<b>ANTI-TIP WHEELS</b>	N/A
<b>MAXIMUM SPEED <sup>4</sup></b>	Up to 6 mph (Optional: 8 mph)
<b>GROUND CLEARANCE</b>	3.25" (frame base)
<b>TURNING RADIUS <sup>5</sup></b>	35.75"
<b>OVERALL LENGTH</b>	36" without foot riggings
<b>BASE WIDTH</b>	25"
<b>SEATING SIZES</b>	
<b>Synergy® / Static Seat</b>	<b>Width:</b> 12-22" <b>Depth:</b> 12-22"
<b>TRU-Balance® 3 Power Tilt</b>	<b>Width:</b> 12-22" <sup>3</sup> <b>Depth:</b> 12-22" <sup>3</sup>
<b>TRU-Balance® 3 Lift &amp; Tilt</b>	<b>Width:</b> 12-22" <sup>3</sup> <b>Depth:</b> 12-22" <sup>3</sup>
<b>TRU-Balance® Flex Seat</b>	<b>Width:</b> 15-22" <b>Depth:</b> 15-20"
<b>SEAT-TO-FLOOR HEIGHTS</b>	
<b>Synergy®/Static Seat <sup>8</sup></b>	Static: 16.5 - 18.1" Synergy: 15.875 - 17.875"
<b>TRU-Balance® 3 Power Tilt <sup>8</sup></b>	18-19.6"
<b>TRU-Balance® 3 Lift &amp; Tilt <sup>8</sup></b>	18.5-20.1"
<b>TRU-Balance® Flex Seat</b>	16.5-18.1"
<b>MANUAL TILT/RECLINE <sup>10</sup></b>	Yes
<b>BATTERY SIZE <sup>6</sup></b>	Group 24
<b>BATTERY WEIGHT <sup>14</sup></b>	52 lbs.
<b>AVAILABLE ELECTRONICS <sup>12</sup></b>	100A Q-Logic 2 EX 75A Q-Logic NE+ 75A Q-Logic NE
<b>BATTERY CHARGER</b>	8A Off-board
<b>MOTOR PACKAGES</b>	4-Pole 6 mph <sup>1</sup> or 8 mph <sup>2</sup>
<b>WEIGHT CAPACITY <sup>7</sup></b>	300 lbs.
<b>BASE WEIGHT</b>	129 lbs.
<b>MEDICARE CODE COVERAGE</b>	<b>Group 3:</b> Standard (K0848/K0849), Single Power (K0856), Multiple Power (K0861)

- 1) Standard
- 2) Optional
- 3) Requires HD Suspension
- 4) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 5) Without front riggings
- 6) Battery capacity may change when adding power positioning systems.
- 7) Weight capacity is dictated by the selection of seating.
- 8) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 9) 250 lbs. weight capacity with TRU-Balance® Systems
- 10) The Synergy® Manual Tilt weight capacity is 75 lbs. and maximum dimensions are 14" W x 14" D. Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.
- 11) Requires Group 34 batteries
- 12) All Quantum Rehab® units can be ordered with specialty controls
- 13) Required with high-speed motors
- 14) Battery weight may vary +/- 3 lbs. based on manufacturer.

# Features

- Compact size for outstanding maneuverability in tight spaces
- Maximum speed up to 8 mph<sup>2</sup>
- High-performance, high-efficiency motor package
- Accepts a wide range of seating and positioning systems
- Sliding battery tray for easy access
- Two suspension packages

# Options

- TRU-Balance® 3 Power Tilt
- TRU-Balance® 3 Power Recline
- TRU-Balance® 3 (10") Lift
- Manual and Power Elevating Leg Rests
- TRU-Comfort 2 Seating
- Q-Logic 2 Controller
- Full Lighting Package
- Centermount Adjustable Foot Platform
- Accu-Trac Advanced Tracking Technology
- Power Articulating Foot Platform

# Colors



Pretty in Pink



Green Machine



Deep Purple



Evolution Orange



Electric Blue



Red Corvette



Grey Street



Midnight Blue



Night Train



Back in Black